

PULPOGEL®

Ferric sulfate gel for vital pulpotomy of primary teeth and bleeding control

EFFECTIVE AND SAFE PRIMARY TEETH TREATMENT



NEW PRODUCT



2,4 g

C € 2282


PULPOGEL[®] application in vital pulpotomy of primary teeth



Make local anesthesia



After cavity preparation remove roof of the pulp chamber



Remove the pulp at the height of the root canal orifices



Rub preparation into bleeding pulp for 15 sec., if necessary repeat procedure



After bleeding control rinse by air and water-flow, and dry



Apply Zinc oxide-Eugenol (ZOE) cement and final restoration

Intended use

- · Preparation on the basis of ferric sulfate is intended for pulpotomy for vital amputation in primary teeth.
- · Preparation is useful for minor bleeding control that is a result of dental treatment.

Technical and Medical Data

Composition: ferric sulfate - 20%, gelling agent, purified water.

Primary package: syringe with protective cap, containing 2,4 g of gel-type formulation product.

Mechanism of action

Gel inducts haemostasis: during contact with blood it forms protein-sulfate complexes, which are closing blood vessels. It's advantage is simplicity and quickness of application: rub it into the pulp horns. Next, in case of vital pulpotomy, after removal from the chamber (by air and water-flow) apply ZOF and final restoration.

Preparation may be used as a haemostatic agent for capillary hemorrhage during intraoral procedures. Strongly rubbed in bleeding tissue, it forms clots of coagulated blood located in orifice of capillary. Rubbing also cuts off part of a clot that is standing out of the vessel's level and forms a cap permanently closing the vessel. Tissues should be cleaned by a strong flush of air and water-flow, which removes clots, remains of gel and verifies effectiveness of haemostasis.

